



# Sectra Viewer Shortcut Keys

## Overview

|                |                              |
|----------------|------------------------------|
| <b>Tab</b>     | Show/hide image overview     |
| <b>i</b>       | Show/hide slide details view |
| <b>l</b>       | Show/hide annotation list    |
| <b>h</b>       | Show/hide annotations        |
| <b>r</b>       | Show/hide review tracking    |
| <b>Shift r</b> | Toggle reviewed state        |
| <b>f</b>       | Toggle Full Screen           |
| <b>?</b>       | Show shortcut keys help      |
| <b>Esc</b>     | Cancel/close                 |

## Annotations

|              |                                |
|--------------|--------------------------------|
| <b>p</b>     | Arrow annotation               |
| <b>m</b>     | Distance measurement           |
| <b>a</b>     | Area annotation                |
| <b>s</b>     | Square annotation              |
| <b>c</b>     | Circle annotation              |
| <b>Enter</b> | Name the selected annotation   |
| <b>Del</b>   | Delete the selected annotation |

## Image Navigation

|                            |  |
|----------------------------|--|
| <b>Shift Drag</b>          | Rotate image                             |
| <b>Ctrl Scroll</b>         | Adjust image focus                       |
| <b>Click Mousewheel</b>    | Mouse lock mode                          |
| <b>Page Up / Page Down</b> | Zoom in/out stepwise                     |
| <b>Num 4</b>               | Pan a whole view left                    |
| <b>Num 6</b>               | Pan a whole view right                   |
| <b>Num 8</b>               | Pan a whole view up                      |
| <b>Num 2</b>               | Pan a whole view down                    |
| <b>Home</b>                | Reset view                               |
| <b>◀ / ▶</b>               | Go to previous/next image                |
| <b>Shift ▶ / ▶</b>         | Go to first image of previous/next block |
| <b>Num 9 / Num 3</b>       | Go to previous/next annotation           |
| <b>Num + / Num -</b>       | Activate/deactivate sync of images       |

